

SAM SPLINT Finger Splint, Pack of 10



SAM SPLINT Finger Splint for Rapid Immobilization

The **SAM SPLINT finger splint** is a compact immobilization device designed to quickly, easily, and effectively immobilize an injured finger. It is particularly well-suited for first aid situations, emergency response kits, medical, sports, and professional settings, as well as for initial on-site care.

Based on the design of classic SAM SPLINT splints, it combines **flexible shaping** with **structured support**. Once molded, the splint stabilizes the finger to limit movement and provide initial support until medical advice or treatment is available.

A compact size for finger injuries

With its specific dimensions of **9.4 x 4.6 cm**, this version is designed for **finger immobilization**. More compact than a long splint, it is easy to transport and can be stored in a first-aid kit, medical bag, or emergency kit without taking up too much space.

Packaged in **sets of 10 splints**, this ensures a supply of ready-to-use units, particularly for professionals, first-aid organizations, sports facilities, local governments, or settings requiring large quantities of first-aid supplies.

A lightweight, moldable, and reusable splint

The SAM SPLINT finger splint is **lightweight, compact, and easily moldable**. It can be folded or shaped as needed to fit the area requiring immobilization. It is easily secured using medical adhesive tape, a bandage, or a suitable support system.

Its **reusable** and **water-resistant** design makes it a practical solution for a variety of emergency response environments. It is also **radiolucent**, allowing X-ray examinations to be performed without significant interference from the splint.

When should you use a SAM SPLINT finger splint?

This splint is intended for **the temporary immobilization of a finger** in the event of an injury, sprain, or suspected fracture. It serves as a first-aid solution to limit movement while awaiting medical evaluation.

It can be used by **first responders**, medical teams, sports facilities, businesses, local governments, or any organization looking to supplement its emergency kit with a simple and reliable solution for finger injuries.

Is this splint suitable for a first-aid kit?

Yes, its compact size and packaging in **packs of 10** make it ideal for first-aid kits, emergency response bags, sports kits, or first-aid stations.

Can the splint be reused?

Yes, this model is designed to **be reusable**. However, you should check its condition after use and follow the appropriate hygiene guidelines for the situation.

How do you secure the splint in place?

The splint can be secured using **medical adhesive tape**, a bandage, or a suitable fastening device. The choice of fastening depends on the situation, the desired level of comfort, and the advice of the healthcare professional or first responder on the scene.